Test Software P/N 4102483

**Tested Part Number** 1054096

Test Software Name Trilogy FSA (PM TEST)

PHILIPS RESPIRONICS

Software Revision 24.0.0.0

Prepared By P. Pascal

## Test Report

Step	Description	Units	Limits	Results	Pass	Fail
1	0010.0003 Write Charger Limiter Table (Limit=75%) [19s]	N/A	Pass	TRUE		<del></del>
2	0010.0020 Check Leak (1) 63 @ 25,0cmH2O: Control Flow Sen:	SLPM	-5 to 5	1,16	<b>√</b>	
3	0010.0021 Check Leak (1) @ 25,0cmH20: Test Setup [4m 3s]	SLPM	-5 to 5	0,00	<b>√</b>	
4	0010.0070 Check HW Revision [4m 3s]	N/A	EQ 0	0,00	<b>√</b>	
5	0010.0080 Check SW Revision (DSP rev. 13) [4m 3s]	N/A	EQ 14.2.05	14.2.05	<b>√</b>	
6	0010.0110 Check Clock Settings [4m 3s]	s	LE 300	12	<b>√</b>	
7	0010.0120 Int.Batt.Cap.@(T=25C,SH=96%,CC=7,CF=0,ME=1) [4m	8	15 to 85	46,00	<b>√</b>	
8	0010.0121 Det.Batt.Cap.@(T=25C,SH=86%,CC=22,CF=0,ME=1) (1	8	15 to 85	60,00	<b>√</b>	
9	0010.0131 Bet.Batt.Cap.w(1-23c,Sh-80%,CC-22,CF-0,MB-1) (1	mV	2450 to 2550	2502	<b>√</b>	
10	0010.0140 Motor Temperature [4m 3s]	C Deg	15 to 82	31	<b>√</b>	
11	0010.0150 Check CPLD SW Revision [4m 3s]			12	<b>√</b>	
		N/A	EQ 12		<b>√</b>	
12	0010.0160 Boot Monitor SW Revision [4m 3s]  0010.0170 Int. Batt. S/N [4m 3s]	N/A N/A	EQ 4.0 0022DFCB	4.0 0022DFCB	<b>√</b>	
14	0010.0170 Int. Batt. S/N [4m 3s]		0022DFCB	0022DFCB	<b>√</b>	
	0020.1030 Write Image Table [4m 14s]	N/A			<b>√</b>	
15		N/A	Pass	TRUE	<b>√</b>	
16	0030.0010 Sensor Board Table Active [5m 32s]	N/A	Pass	TRUE	<b>√</b>	
17	0030.0020 Device Table Active [5m 32s]	N/A	Pass	TRUE	<b>√</b>	
18	0030.0030 Proximal Pressure Table Active [5m 32s]	N/A	Pass	TRUE	✓	
19	0030.0040 Charger Settings Table Active [5m 32s]	N/A	Pass	TRUE	<b>√</b>	
20	0030.0050 Image Table Active (Trilogy100_ImageCalTable.bi	N/A	Pass	TRUE	<b>√</b>	
21	0030.0080 Device Name [5m 35s]	N/A		ogy 100, Internati	<b>√</b>	
22	0030.0090 Device Model [5m 35s]	N/A	EQ 1054096	1054096	✓	
23	0030.0100 Device S/N [5m 35s]	N/A	EQ TV120092328	TV120092328	✓	
24	0030.0110 Product ID (Trilogy 100 Ventilator) [5m 35s]	N/A	EQ 2C	2C	✓	
25	0030.0170 Pos Flow Verify: dP2 at 189,3 (Setpoint 190) [1	SLPM	178,3 to 200,4	195,3	✓	
26	0030.0180 Pos Flow Verify: dP2 at 166,5 (Setpoint 165) [1	SLPM	156,3 to 176,7	173,2	✓	
27	0030.0190 Pos Flow Verify: dP2 at 142,8 (Setpoint 140) [1	SLPM	133,6 to 152,0	147,4	✓	
28	0030.0200 Pos Flow Verify: dP2 at 132,2 (Setpoint 130) [1	SLPM	123,4 to 141,0	136,9	✓	
29	0030.0210 Pos Flow Verify: dP2 at 120,7 (Setpoint 120) [1	SLPM	112,3 to 129,0	126,8	✓	
30	0030.0220 Pos Flow Verify: dP2 at 107,0 (Setpoint 110) [1	SLPM	99,2 to 114,8	113,5	✓	
31	0030.0230 Pos Flow Verify: dP2 at 98,2 (Setpoint 100) [10	SLPM	90,7 to 105,6	103,4	✓	
32	0030.0240 Pos Flow Verify: dP2 at 87,2 (Setpoint 90) [10m	SLPM	80,2 to 94,2	92,9	✓	
33	0030.0250 Pos Flow Verify: dP2 at 78,0 (Setpoint 80) [10m	SLPM	71,4 to 84,6	82,9	✓	
34	0030.0260 Pos Flow Verify: dP2 at 67,2 (Setpoint 70) [10m	SLPM	61,0 to 73,4	71,4	✓	
35	0030.0270 Pos Flow Verify: dP2 at 58,0 (Setpoint 60) [10m	SLPM	52,2 to 63,8	61,1	✓	

**Test Started On** 04/01/24 07:39:12

Elapsed Test Time 48m 58s

Serial Number  $\underline{\text{TV120092328}}$ 

Status PASS

**Page** 1 **of** 6